



Date: 14.10.2018

Created by: Fabian Kleene

This document is the legal and intellectual property of MK-dent GmbH and may only be used as permitted. Information are solely extracted from the manufacturers catalogue and/or website and the date mentioned above. This document may not reflect the current state as changes could have occurred.

Product name	Model	Power	Sprayholes	Noise level	Glassrod	Carbide insert	Holding force	Integrated Pressure regulation	Material	Warranty
Premium class standard-head size										
MK-dent	Prime Line HP21KL	25	4	<55	Yes	Yes	35N	Yes	Brass	2 Years Unconditional
Kavo	Mastertorque M9000L	23	4	57	Yes	Yes	30N	Yes	Steel	2 Years
Sirona	T1 Control	22	4	60	Yes	No	26N	No	Titanium	2 Years
Sirona	T1 Boost	23	4	60	Yes	No	26N	No	Titanium	2 Years
NSK	Ti-Max Z900L	26	4	62	Yes	No	23N	No	Titanium	2 Years
W&H	Synea Vision TK-100	26	5	57	Yes	No	30N	No	Brass	2 Years
W&H	Synea Vision TK-98	24	5	59	Yes	No	30N	No	Brass	2 Years
B.A. International	Ultimate Power+	20	3	60	Yes	No	26N	No	Titanium	2 Years
Premium class small-head size										
MK-dent	Prime Line HP22KL	18	3	60	Yes	Yes	35N	Yes	Brass	2 Years Unconditional
KaVo	Mastertorque M8700L	19	4	59	Yes	Yes	32N	Yes	Brass	2 Years
Sirona	T1 Mini	20	4	60	Yes	No	26N	No	Titanium	2 Years
W&H	Synea Vision TK-94	14	5	59	Yes	No	30N	No	Brass	2 Years
W&H	Synea Vision TK-97	21	5	59	Yes	No	30N	No	Brass	2 Years
B.A. International	Ultimate Power+	14	3	60	Yes	No	26N	No	Titanium	2 Years
Middle class standard head size										
MK-dent	Classic Line HC21KL	22	4	60	Yes	Yes	35N	Yes	Brass	1,5 Years unconditional 2 Years regular
MK-dent	Classic Line HC21K	22	4	60	No	Yes	35N	Yes	Brass	1,5 Years unconditional 2 Years regular
MK-dent	Classic Line HC20KL	28	4	60	Yes	Yes	35N	Yes	Brass	1,5 Years unconditional 2 Years regular
MK-dent	Classic Line HC20K	28	4	60	No	Yes	35N	Yes	Brass	1,5 Years unconditional 2 Years regular
NSK	X-Max X700L	21	4	65	Yes	No	>25N	No	Titanium	2 Years
NSK	X-Max X700	21	4	65	No	No	>25N	No	Titanium	2 Years
B.A. International	Ultimate	20	3	60	Yes	No	k.A.	No	Brass	1 Year
Middle class small head size										
MK-dent	Classic Line HC22KL	18	3	65	Yes	Yes	35N	Yes	Brass	1,5 Years unconditional 2 Years regular
MK-dent	Classic Line HC22K	18	3	65	No	Yes	35N	Yes	Brass	1,5 Years unconditional 2 Years regular
NSK	X-Max X600L	17	4	65	Yes	No	>25N	No	Titanium	2 Years
NSK	X-Max X600	17	4	65	No	No	>25N	No	Titanium	2 Years
B.A. International	Ultimate	14	3	60	Yes	No	k.A.	No	Brass	1 Year
Entrance class standard head size										
MK-dent	Eco Line HE21KL	24	4	60	Yes	No	35N	No	Brass	1 Year unconditional 2 Years regular
MK-dent	Eco Line HE21K	24	4	60	No	No	35N	No	Brass	1 Year unconditional 2 Years regular
KaVo	Expert E680/E679	18	4	62	Yes	Yes	30N	Yes	Brass	1,5 Years unconditional 2 Years regular
Sirona	T2 Control	22	4	58	Yes	No	k.A.	No	Brass	2 Years
Sirona	T2 Boost	23	4	58	Yes	No	k.A.	No	Brass	2 Years
NSK	S-Max M600L	17	4	62	Yes	No	k.A.	No	Stahl	1,5 Years
NSK	S-Max M600	17	4	62	No	No	k.A.	No	Stahl	1,5 Years

Product name	Model	Power	Sprayholes	Noise level	Glassrod	Carbide insert	Holding force	Integrated Pressure regulation	Material	Warranty
W&H	Synea Fusion TG-98	21	4	59	Yes	No	k.A.	No	Brass	1 Year
B.A. International	Optima	17	3	62	Yes	No	k.A.	No	Brass	1 Year
B.A. International	Optima	17	3	62	No	No	k.A.	No	Brass	1 Year
B.A. International	Optima 525	20	4	65	Yes	Yes	k.A.	No	Brass	1 Year
B.A. International	Optima 525	20	4	65	No	Yes	k.A.	No	Brass	1 Year
Entrance class small head size										
MK-dent	Eco Line HE22KL	16	3	62	Yes	No	35N	No	Brass	1 Year
MK-dent	Eco Line HE22K	16	3	62	No	No	35N	No	Brass	1 Year
KaVo	Expert E677/E675	15	3	62	Yes	Yes	30N	Yes	Brass	1,5 Years
NSK	S-Max M500L	16	4	65	Yes	No	k.A.	No	Stahl	1,5 Years
NSK	S-Max M500	16	4	65	No	No	k.A.	No	Stahl	1,5 Years
Sirona	T2 Mini	20	3	60	Yes	No	k.A.	No	Brass	2 Years
W&H	Synea Fusion TG-97	18	4	59	Yes	No	k.A.	No	Brass	1 Year
B.A. International	Optima 525	18	4	65	Yes	Yes	k.A.	No	Brass	1 Year
Star	430	13	1	69	Yes	No	k.A.	No	Brass	0,5 Years
Midwest	Stylus	11	4	60	Yes	No		No	Aluminium	1 Year
Basic class standard head size - Fixed End 4-hole										
MK-dent	Eco Line HB21KL	21	4	60	Yes	No	35N	No	Brass	1 Year
MK-dent	Eco Line HB21K	21	4	60	No	No	35N	No	Brass	1 Year
MK-dent	Eco Line HE12	18	1	65	No	No	35N	No	Brass	1 Year
MK-dent	Eco Line HB12	18	1	65	No	No	35N	No	Brass	1 Year
MK-dent	Eco Line HE17	28	4	65	No	No	35N	No	Brass	1 Year
MK-dent	Eco Line HB17	18	4	65	No	No	35N	No	Brass	1 Year
KaVo	SMARTorque LUX S619L	16	4	65	Yes	No	30N	No	Brass	1 Year
KaVo S	SMARTorque S619C	16	4	65	No	No	30N	No	Brass	1 Year
KaVo	SMARTorque S609C	15	4	65	No	No	k.A.	No	Brass	1 Year
Sirona	T3 Boost	23	4	62	No	No	k.A.	No	Brass	1 Year
Sirona	T3 Racer	30	4	65	No	No	k.A.	No	Brass	1 Year
Sirona	T4 Racer	17	4	65	No	No	k.A.	No	Brass	1 Year
Midwest	Quiet-air-L Power Lever	10	1	75	Yes	No		No	Aluminium	0,5 Years
NSK	Pana Max Plus	11	4	68	No	No	k.A.	No	Stahl	1 Year
NSK	Pana Max2	11	1	65	No	No	k.A.	No	Stahl	1 Year
Basic class small head size - Fixed End 4-hole										
MK-dent	Eco Line HB22KL	16	3	62	Yes	No	35N	No	Brass	1 Year
MK-dent	Eco Line HB22K	16	3	62	No	No	35N	No	Brass	1 Year
MK-dent	Eco Line HE12	18	1	65	No	No	35N	No	Brass	1 Year
MK-dent	Eco Line HB14	18	1	65	No	No	35N	No	Brass	1 Year
MK-dent	Eco Line HE18	18	4	65	No	No	35N	No	Brass	1 Year
MK-dent	Eco Line HB18	18	4	65	No	No	35N	No	Brass	1 Year
KaVo	SMARTorque mini S615L	13	1	65	Yes	No	30N	No	Brass	1 Year
KaVo	SMARTorque S605C	13	1	65	No	No	k.A.	No	Brass	1 Year
NSK	Pana Max Plus	11	4	68	No	No	k.A.	No	Stahl	1 Year
Midwest	Tradition	13	1	77	Yes	No		No	Aluminium	0,5 Years
Midwest	Quiet-air-L Power Lever	10	1	75	Yes	No		No	Aluminium	0,5 Years
Entrance class standard head size - Fixed End 2-hole										
MK-dent	Eco Line HE11	18	1	65	No	No	35N	No	Brass	1 Year

